


Ultimaster Nagomi[™]


Sirolimus eluting coronary stent system

READY FOR ACTION.



Optimized overexpansion capability
up to $\varnothing 6.25 \text{ mm}^3$

NEW **10-crown design*** for a better
expansion in large vessels^{1,2,3}



NEW **Balloon with a nominal pressure
of 11 atm** for a better apposition^{4,5}

www.terumo-europe.com

* for $\varnothing 3.50 - 4.50 \text{ mm}$ stents

- 1 Technical documentation, PS-0084
- 2 Product Standard, PS-1084
- 3 Ultimaster Nagomi[™] Instructions for Use
- 4 Internal test report, LMF-0000041-RAFR-001-1070
- 5 Internal test report, LMF-0000041-RAFR-001-1159

Ultimaster Nagomi[™] is not available for sale in all countries.
Please contact your Terumo local sales representative for more information.